

# KR-1

Auto Kerato-Refractometer



 **TOPCON Healthcare**

SEEING EYE HEALTH DIFFERENTLY

# Fully automated Easy-to-use Kerato-Refractometer




## Overview



**Automated  
Kerato-Refractometer**



**Rotary Prism  
Technology**



**Full 360° rotating monitor**  
allows operator to  
maintain distance from patient



Small footprint.  
**Space Saving**



**Patient Friendly**



**Easy to use**

## KR-1

The Topcon KR-1 is a modern, automated, easy-to-use kerato-refractometer. Controlled from the color touchscreen and incorporating Topcon's Rotary Prism technology.



# Automatic Workflow



## Simple Operation

Simply touch the center of the pupil on the touchscreen control panel, and the KR-1 automatically captures right then left eye information. The KR-1 provides an innovative operator experience.



Touch the pupil image on the screen for auto-alignment



Automatic capture and R/L movement



Data displayed with option to transfer and/or print

# KR-1 Features

## Easy-to-Use Color Touchscreen Control Panel

All operating procedures can be performed with the touchscreen control panel. The optimized user interface provides an intuitive operator experience. The KR-1 rapidly captures data from the right eye before moving automatically to capture the left eye. An easy-to-use manual capture option is also available, featuring manual focus coupled with auto-alignment.



# Flexible and Space Saving Layout

The control panel can be rotated in any direction. The compact device enables the operator to support the patient easily during measurement, even if the patient's eyelid needs to be lifted. These unique aspects will contribute to space saving and flexible layout in your eye examination room.



Operation from  
the front  
(Seating)



Operation from  
the front  
(Standing)



Operation from  
the side




Operation from  
the patient's side

# KR-1 Features

## Compact Ergonomic Design

The KR-1's ergonomic design provides a more comfortable position for the patient with an inclined optical head (5°). The compact instrument enables easier access for patients.




## Peripheral KRT

The KR-1 allows the operator to capture 8 peripheral corneal curvature measurements in addition to the central data.

## Topcon's Rotary Prism Measuring System

With Topcon's Rotary Prism Technology, the KR-1 provides stable measurement. The unique, eccentric rotation of the measurement ring, is designed to reduce the artifact from the fundus.

## System Chart



\*Please contact Topcon subsidiaries in your region for the detailed system configurations.

# Specification of KR-1

<b>Refractive Power Measurement</b>	Spherical refractive power: -25D to +22D (0.12D/0.25D steps)*
	Cylindrical refractive power: 0D to ±10D (0.12D/0.25D steps)*
	Astigmatic axis angle: 0° to 180° (1°/ 5°steps)
	Minimal measurable pupil diameter: φ 2 mm
<b>Corneal Curvature Measurement</b>	Corneal curvature radius: 5.00mm to 10.00 mm (0.01mm step)
	Corneal refractive power: 67.50D to 33.75D (0.12D/0.25D steps) (keratometric corneal refractive index = 1.3375)
	Corneal astigmatic axis angle: 0° to 180°(1°/ 5° steps)
	Corneal astigmatic refractive power: 0D to ±10D (0.12D/0.25D steps)
<b>PD Measurement Range</b>	20mm to 85mm (1mm step)
<b>Data Transport Terminal</b>	USB (Import), RS-232C (Export), LAN (Export)
<b>Dimensions</b>	286-326mm (W) × 445-526mm (D) × 466-615 mm (H)
<b>Weight</b>	19kg
<b>Power Supply</b>	100-240V AC, 50-60Hz, 75VA

\* -25D ≤ spherical refractive power + cylindrical refractive power and spherical refractive power + cylindrical refractive power ≤ +22D

Not all products, services or offers are approved or offered in every market, and products vary from one country to another. Contact your local distributor for country-specific information and availability.

**IMPORTANT** In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation.

## TOPCON CORPORATION

75-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8580, JAPAN. Phone: +81-(0)3-3558-2522/2502 Fax: +81-(0)3-3965-6898 | topconhealthcare.jp

### TOPCON SINGAPORE MEDICAL PTE. LTD.

100G Pasir Panjang Road, #05-05, Interlocal Centre, SINGAPORE 118523  
Phone: +65-68720606 Fax: +65-67736150  
E-mail: med.sales.sg@topcon.com  
topconhealthcare.sg

### TOPCON INSTRUMENTS ( MALAYSIA ) SDN. BHD.

No.6, Jalan Pensyarah U1/28, Hicom Glenmarie Industrial Park, 40150 Shah Alam, Selangor, MALAYSIA  
Phone:+60-(0)3-50223688 Fax: +60-(0)3-50313968  
E-mail: mys\_medical\_sales@topcon.com  
topconhealthcare.my

### TOPCON INSTRUMENTS ( THAILAND ) CO., LTD.

77/162 Sinnsathorn Tower, 37th Floor, Krungthonburi Rd., Klongtongsai, Klongsarn, Bangkok 10600, THAILAND  
Phone: +66(0)2-440-1152-7 Fax: +66-(0)2-440-1158  
E-mail: tha\_medical@topcon.com  
eyecare.topcon.co.th

### MEHRA EYETECH PRIVATE LIMITED

801 B Wing, Lotus Corporate Park, Graham Firth Steel Compound Goregaon (East) Mumbai 400063 Maharashtra, INDIA  
Phone: +91-22-61285455  
E-mail: sales@mehraeyetech.in  
topconhealthcare.in

### TOPCON (BEIJING) MEDICAL TECHNOLOGY CO., LTD.

Room 2808, Tower C, JinChangAn Building, No.82, Middle Section of East 4th Ring Road, Chaoyang District, Beijing 100124, P.R. CHINA

